



CALIBRE PCM

Instrumento de sobremesa adecuado para determinar la longitud, ancho y diagonal de baldosas cerámicas o productos similares.

Se puede suministrar solo la parte superior para colocar sobre un banco o sobre un marco específico (se debe pedir por separado)

Especificaciones técnicas

- Display LCD con resolución 0,1mm
- Precisión +/-0,1mm
- Sistema de medición magnético

Dotación:

- 2 pilas tipo AAA
- Estructura de acero pintado
- Plano de soporte de baldosa en acero inoxidable

PCM CALIBER

Bench-top instrument suitable for determining the length, width, diagonal of ceramic tiles, or similar products.

Can be supplied only the top to be placed on a bench or on a dedicated frame (to be ordered separately)

TECHNICAL SPECIFICATIONS

- LCD Display with resolution of: 0,1 mm
- Precision $\pm 0,1$ mm
- Magnetic measuring system

EQUIPMENT

- Nr 2 AAA-Type batteries
- Painted steel structure
- Support plan of tile in stainless steel

Código / Code	Modelo / Model	Medida palpador	Dimensiones externas	Peso / Weight
GT2179	PCM650	650mm	890x750x170mm	60Kg
GT2182	PCM W-650	650mm	890x750x170mm	60Kg
GT2180	PCM1250	1250 mm	1490x750x170mm	110Kg
GT2183	PCM W-1250	1250mm	1490x750x170mm	110Kg
GT2181	PCM1500	1500mm	1740x750x170mm	160Kg
GT2184	PCM W-1500	1500mm	1740x750x170mm	160Kg